

U.S. ARMY RESEARCH, DEVELOPMENT AND ENGINEERING COMMAND

ECBC-TR-662

# EXPLORATION OF THE USE OF NUCLEAR MAGNETIC RESONANCE FOR THE STUDY OF RICIN TOXICITY IN CELLS

Vicky L. H. Bevilacqua Janna S. Madren-Whalley Lisa M. Reilly Saumil S. Shah Jeffrey S. Rice

RESEARCH AND TECHNOLOGY DIRECTORATE

April 2009

20090526375

Approved for public release; distribution is unlimited.



Disclaimer
The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorizing documents.

data needed, and completing this burden to Department of 4302. Respondents should be	and reviewing this collection of it Defense, Washington Headquart e aware that notwithstanding any	nformation. Send comments rega ers Services, Directorate for Infon	rding this burden estimate or any mation Operations and Reports ( shall be subject to any penalty f	other aspect of this co 0704-0188), 1215 Jeffe	hing existing data sources, gathering and maintaining the illection of information, including suggestions for reducing erson Davis Highway, Suite 1204, Arlington, VA 22202- n a collection of information if it does not display a currently	
1. REPORT DATE (DI XX-04-2009	D-MM-YYYY)	2. REPORT TYPE Final	E33.		PATES COVERED (From - To) ar 2008 - Oct 2008	
4. TITLE AND SUBTIT		Magnetic Resonar	nce for the Study o		CONTRACT NUMBER	
Exploration of the Use of Nuclear Magnetic Resonan Toxicity in Cells			ioo ioi iiio olaay oi riioiii		b. GRANT NUMBER	
				5c.	PROGRAM ELEMENT NUMBER	
6. AUTHOR(S Bevilacqua, Vicky L. H.; Madren-Whalley, Janna S.; Saumil S.; and Rice, Jeffrey S.			Railly Lisa M · Sh		PROJECT NUMBER 61110191A00	
			rtelliy, Lisa W., Orlan,		TASK NUMBER	
				5f. \	WORK UNIT NUMBER	
	GANIZATION NAME(S) TN: AMSRD-ECB	AND ADDRESS(ES) -RT-DD, APG, ME	21010-5424	N	PERFORMING ORGANIZATION REPORT JUMBER ECBC-TR-662	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(			(ES)	10.	SPONSOR/MONITOR'S ACRONYM(S)	
					SPONSOR/MONITOR'S REPORT NUMBER(S)	
	AVAILABILITY STATEM blic release; distri	ENT bution is unlimited.				
13. SUPPLEMENTAR	YNOTES					
Laboratory Indepleased solely on effects. Comple of statistical methas required moon an explorator BALB/c 3T3 Multiple cell growth intool. We have for same preparation.	pendent Research cell survivability. ex cellular metabol hods such as che del-building with lay study to determine Fibroblasts (Conedia alone. This bund that NMR-Poen that had cells grence that this process.	n project "Cell Toxi However, addition ite mixtures can be mometrics analysi C <sub>50</sub> /EC <sub>50</sub> values of ine the feasibility of CCL-163 American work involves prin CA readily distingurown in them from	city by NMR." Tra al information is de e analyzed by nuc s (CA). However, otained from conve of using NMR-CA i Type Culture Coll acipal component a ishes between mu those that did not	esirable for user magnetic the use of Nentional cytos ndependenticetion [ATCO analysis (PC altiple-compo have cells greational assertions of the collection for the cells greation for the cells greating analysis of the	Biological Center In-House ays lead to IC <sub>50</sub> /EC <sub>50</sub> values understanding a toxin's total c resonance (NMR) with the aid IMR-CA for cell toxicity to date toxicity assays. Here we report y for the study of ricin toxicity in C], Manassas, VA) by analysis of A) as a chemometrics analysis ment media samples from the rown in them. We also have a media from cells dosed with	
15. SUBJECT TERMS 3T3 Cells <i>Ricinus commur</i>	Cell Tox	cicity Nuclear	Magnetic Resona	ance	NMR Ricin	
16. SECURITY CLAS	SIFICATION OF:		17. LIMITATION OF ABSTRACT	18. NUMBER OF PAGES	19a. NAME OF RESPONSIBLE PERSON Sandra J. Johnson	
a. REPORT U	b. ABSTRACT U	c. THIS PAGE U	UL	14	19b. TELEPHONE NUMBER (include area code) (410) 436-2914	
					Standard Form 298 (Rev. 8-98) Prescribed by ANSI Std. Z39.18	

REPORT DOCUMENTATION PAGE

Form Approved OMB No. 0704-0188 Blank

#### **PREFACE**

The work described in this report was authorized under Project No. 61110191a00, In-House Laboratory Independent Research. This work was started in March 2008 and completed in October 2008.

The use of either trade or manufacturers' names in this report does not constitute an official endorsement of any commercial products. This report may not be cited for purposes of advertisement.

This report has been approved for public release. Registered users should request additional copies from the Defense Technical Information Center; unregistered users should direct such requests to the National Technical Information Service.

#### **Acknowledgments**

The authors would like to acknowledge the In-House Laboratory Independent Research program for funding. The authors would also like to acknowledge Dr. H. Dupont Durst, Jude Height, and Alan Zulich for administrative support.

Blank

## CONTENTS

1.	INTRODUCTION	
2.	METHODS	
2.1 2.2 2.3 2.4 2.5 2.6	BALB/c 3T3 Murine Fibroblasts Ricin Preparation Ricin Exposures NMR Sample Preparation NMR Experiments Chemometrics Analysis	8 8
3.	RESULTS AND DISCUSSION	9
3.1 3.2	NMR ExperimentsPCA Analysis	9
4	CONCLUSIONS	10

# **FIGURES**

ı	Sample that Contains Media on which Cells Were Grown	
2	PC2 versus PC1 Scores Plot for PCA Involving 76 Spectra	,
3	PC3 versus PC4 Scores Plot for PCA Involving 76 Spectra	)
4	PC2 versus PC1 Scores Plot for PCA Involving Batch 1 Spectra	•
5	PC2 versus PC1 Scores Plot for PCA Involving Batch 4 Spectra	
	TABLE	
	Summary of Samples analyzed by NMR9	)

# EXPLORATION OF THE USE OF NUCLEAR MAGNETIC RESONANCE FOR THE STUDY OF RICIN TOXICITY IN CELLS

#### 1. INTRODUCTION

Traditional cell toxicity assays lead to an  $IC_{50}$  or  $EC_{50}$  value based on the determination of cell death. A method that provides information in addition to survivability regarding the effect of a toxin on the cell is desirable for understanding the total effects of the toxin. Complex cellular metabolite mixtures can be analyzed by nuclear magnetic resonance (NMR) with the aid of statistical methods such as chemometrics analysis (CA). When combined with CA, NMR may provide information in the form of biomarkers, affording detailed biochemical information not available through typical cell assays, while also allowing for the determination of  $EC_{50}$  or  $IC_{50}$ . An NMR-based CA procedure is also expected to have advantages including: no dyes or labeled compounds; the analysis of dosed cells with little sample preparation; and no additional washes or incubations after the dosing period. Here we report on an exploratory study to determine the feasibility of using NMR-CA independent of conventional cytoxicity assays for the study of ricin toxicity in BALB/c 3T3 Murine Fibroblasts (CCL-163 American Type Culture Collection [ATCC], Manassas, VA) by analysis of the cell growth media alone. This work involves principal component analysis (PCA) as the chemometrics analysis tool.

#### METHODS

Note: Ricin is extremely toxic. Its use is controlled under the Biological Select Agents and Toxins program in the United States by the Centers for Disease Control and Prevention (Atlanta, GA). Handling of ricin should follow strict safety procedures determined in collaboration with the safety office of the research laboratory's organization.

#### 2.1 BALB/c 3T3 Murine Fibroblasts.

The fibroblasts (CCL-163 American Type Culture Collection [ATCC], Manassas, VA) were passaged a minimum of three times after thawing prior to testing. Ninety-six well plates were seeded at  $2.5 \times 10^3$  cells per well and maintained in culture at 37 °C in a humidified atmosphere of 5% CO<sub>2</sub> in air for 24 hr prior to treatment.

#### 2.2 <u>Ricin Preparation</u>.

The *Ricin communis agglutinin* II (ricin) stock solution was prepared by dialyzing ricin (Vector Laboratories, Burlingame, CA) into 10 mM sodium phosphate buffer (pH 7.0, PB) over a period of ~ 24 hr with gentle stirring using three volumes of ~ 600 mL PB. Dialysis was carried out on ice using regenerated cellulose or cellulose ester Dispodialyzers® (Spectrum Laboratories, Rancho Dominguez, CA) with 5000 or 8000 molecular weight cutoff. Ricin stock solutions were stored at 0–4 °C. Stock solution concentration was determined at 25 °C by absorbance of a stock solution dilution using a JASCO Model J-810 Spectropolarimeter (JASCO Analytical Instruments, Easton, MD) equipped with a Peltier thermoelectric temperature control system (PTC-423S). Fifteen absorbance measurements on the solvent solution without any protein were recorded at intervals of 1 s and averaged. Fifteen measurements were then

recorded on the sample and averaged. The blank average was subtracted from the sample average. The concentration was calculated using Beer's Law with  $E^{0.1\%}_{280nm} = 1.4^{1}$ . The ricin was sterile filtered prior to use in the cell cultures.

#### 2.3 Ricin Exposures.

Exposures were performed  $24 \pm 2$  hr after seeding the 96 well plates. Eight different concentrations of ricin were prepared by diluting the ricin in cell culture medium. The diluted ricin was then added to 6 wells (per concentration) in each of 8 plates (Ricin Cells). These plates also included Vehicle Control wells that contained cells to which a 1:1 ratio of growth medium and chemical dilution medium was added (i.e., the wells contained cells, but no ricin). After ricin exposure, the plates were incubated at 37 °C in a humidified atmosphere of 5% CO<sub>2</sub> in air for 48 hr. Four additional 96 well plates (Ricin Blanks) that had not been seeded with cells also received ricin concentrations identical to those used to expose the plates containing cells. These plates also contained Media Blanks (no ricin, no cells). These plates were incubated in the same manner as those with cells.

#### 2.4 <u>NMR Sample Preparation</u>.

Collection of NMR samples was done 48 hr ± 0.5 hr post-exposure. A total of 76 samples were collected. For each Ricin Cells sample, medium from each of the 6 wells of a given ricin concentration from one plate (8 total concentrations per plate) were pooled along with two additional wells (total of 8 wells) from duplicate plates to provide adequate volume for one NMR sample. A total of 6 samples for each ricin concentration were collected (Ricin Cells) along with six samples of Vehicle Control exposed medium. Corresponding samples (two samples per concentration, 8 concentrations) were also collected using the four 96 well plates that were exposed to ricin, but that contained no cells (Ricin Blanks), along with six Media Blank samples. After collection, all samples were filtered using a 0.2 µm sterile syringe filter to remove any possible remaining cell debris from the pooled medium. Samples were quickfrozen on liquid nitrogen and stored at -80 °C. Frozen samples were later thawed to room temperature, concentrated to near-dryness using a Centrivap (Labconco, Kansas City, MO), and stored at -20 °C until NMR analysis. The samples were separated into 4 batches (Batches 1-4) of 19 samples each for drying because the centrivap rotor could not accommodate them all. Each Batch contained samples representing the range of ricin concentrations from the Ricin Cells samples, samples representing at least three of the ricin concentrations from the Ricin Blank samples, at least one Vehicle Control sample, and at least one Media Blank sample.

#### 2.5 NMR Experiments.

On the day of analysis, each sample was dissolved in  $602~\mu L~D_2O$  containing 0.5 mM 2,2-dimethyl-2-silapentane-5-sulfonic acid (DSS) and 0.01 % NaN<sub>3</sub> (chemicals from Sigma-Aldrich, St. Louis, MO). Samples were randomly taken from differing Batches on different days of analysis. One-dimensional proton NMR experiments were carried out on a Varian UNITY Inova 600 MHz NMR spectrometer at a regulated temperature of 25 °C. Experiments used 16K complex points, 1.00 s relaxation delay, 2.37 s acquisition time, and 6913.8 Hz spectral width. Water suppression was carried out with a soft pulse at the water resonance for 2.00 s.

\_

<sup>&</sup>lt;sup>1</sup> Shimoda, T.; Funatsu, G. Binding of Lactose and Galactose to Native and Iodinated Ricin D. *Agric. Biol. Chem.* **1985**, *49*, 2125-2130.

### 2.6 <u>Chemometrics Analysis</u>.

The NMR spectra were initially processed with *NUTS*<sup>2</sup> software. The NMR free induction decays (FIDs) were converted to the frequency domain using fast Fourier transformation. The transformed spectra were individually phased, baseline corrected, and referenced to DSS at 0.00 ppm. The data was reduced by segmenting the spectra from 0.2 to 10 ppm into 0.04 ppm regions, which were integrated. The water resonance region of 4.5 -5.0 ppm was excluded. The resulting 233 values were normalized to allow for comparison among the different spectra. The data sets were exported to *The Unscrambler*<sup>3</sup> software for principal component analysis (PCA) and cluster analysis. Prior to analysis, the data sets were scaled using mean center scaling.

#### RESULTS AND DISCUSSION

#### 3.1 NMR Experiments.

Proton NMR experiments were carried out on 76 samples as summarized in the Table below. A typical 1-dimensional NMR spectrum obtained for cell media on which samples were grown, is shown in Figure 1.

Table. Summary of Samples analyzed by NMR

		# of
Sample Type	Ricin Concentration	Duplicate
	(μg/mL)	Samples
Ricin Cells: media, cells grown, ricin		
Concentration 1	6.0000 x 10 <sup>-03</sup>	6
Concentration 2	4.0816 x 10 <sup>-03</sup>	6
Concentration 3	2.7766 x 10 <sup>-03</sup>	6
Concentration 4	1.8889 x 10 <sup>-03</sup>	6
Concentration 5	1.2849 x 10 <sup>-03</sup>	6
Concentration 6	8.7410 x 10 <sup>-04</sup>	6
Concentration 7	5.9463 x 10 <sup>-04</sup>	6
Concentration 8	4.0451 x 10 <sup>-04</sup>	6
Vehicle Controls: media, cells grown, no ricin	0	6
Ricin Blanks: media, no cells grown, ricin		
Concentration 1	6.0000 x 10 <sup>-03</sup>	2
Concentration 2	4.0816 x 10 <sup>-03</sup>	2
Concentration 3	2.7766 x 10 <sup>-03</sup>	2
Concentration 4	1.8889 x 10 <sup>-03</sup>	2
Concentration 5	1.2849 x 10 <sup>-03</sup>	2
Concentration 6	8.7410 x 10 <sup>-04</sup>	2
Concentration 7	5.9463 x 10 <sup>-04</sup>	2
Concentration 8	4.0451 x 10 <sup>-04</sup>	2
Media Blanks: media, no cells grown, no ricin	0	6

<sup>&</sup>lt;sup>2</sup> NUTS, version 2D; Acorn NMR, Inc. Livermore, CA. 2004.

<sup>3</sup> The Unscrambler, version 9.6; CAMO Software AS. Oslo, Norway, 2006.

#### 3.2 PCA Analysis.

In brief, PCA calculates a set of descriptors, the principal components (PCs) that are linear combinations of the original variables, the parts per million binned regions in this case, with each PC orthogonal to the others. The first PC, PC1, explains the most variance in the data set whereas each successive PC explains the most variance remaining after the previous PC. A plot of PC1 versus PC2 is a two dimensional representation of the complete data set. A graphical separation between groups of data points on a PC scores plots indicates a toxicity-induced variation between the groups of spectra. Differences may be a result of multiple changes in the NMR spectra. The ppm regions responsible for the variation of the data points can be identified using a second graphical tool, the loadings plot.

An initial PCA analysis was carried out using all 76 samples. The analysis was carried out on binned spectra with a bin width of 0.04 ppm and covered the 0.20 – 10 ppm spectral region. The water peak region (4.5 to 5.0 ppm) was excluded from the calculations. Figure 2 shows the PC2 versus PC1 scores plot. No outliers were removed for this first set of calculations. Other than the obvious outliers (red "RC" points, left side of plot), most noticeable is a clear separation based on the Centrivap Batch number. This is accounted for by the fact that the Centrivap developed a vacuum leak that was not initially noticeable, but that increased over time, with the result that the dry-down time successively increased for each Batch dried until the leak was noticed. The separation in the scores plot follows a trend corresponding to the length of dry-down time (Batch 4 = shortest, Batch 3 = longest).

Figure 3 shows the PC4 versus PC3 scores plot from the same analysis (all 76 samples, no outliers removed). Here, separation occurs between the media in which cells were grown and the media in which no cells were grown but which were otherwise incubated and processed in the same manner. Because of the obvious effect of the dry-down time on the analysis noted in the PC2 versus PC1 discussion above. PCA was next carried out on individual batches to see if the separation noted here in the PC4 versus PC3 scores plot would improve. Figure 4 shows the resulting PC2 versus PC1 scores plot for the PCA of Batch 1, this time with outliers removed. The separation is now at the PC2/PC1 level. The lactate chemical shift is the primary contributor to the separation. The concentration of lactate (1.33 ppm NMR doublet, see Figure 1) is higher relative to other metabolites in the media that had cells grown in it compared to the Blanks. The fact that the Blanks with and without ricin do not separate from each other indicates that the separations relate to cell activity rather than to the low levels of ricin that may be present in a sample. To complete the exploratory study, PCA was again carried out on individual Batches, this time with just the Ricin Cells and Vehicle Control samples. Figure 5 shows the resulting PC2 versus PC1 plot for Batch 4. For concentrations having more than one duplicate in the batch, the samples are grouped relative to the ricin-dosing concentration. This result is promising and indicates that there may be metabolic effects that depend on the dosing concentration.

#### 4. CONCLUSIONS

We have found that nuclear magnetic resonance-principal component analysis (NMR-PCA) readily distinguishes between multiple-component media samples from the same preparation which had cells grown in them from those which did not have cells grown in them. In addition, we have shown that NMR-PCA detects variation in the time taken to dry down samples for NMR analysis. Finally, we have preliminary evidence that this procedure will be able to identify differences between media from cells dosed with varying concentrations of ricin.

Additional studies should include a replicate of the study carried out here, but with all samples dried down simultaneously and as quickly as possible. This change in procedure will allow for more statistically relevant PCA from which variation of specific metabolites would be expected to lead to an  $EC_{50}$  determination.




Figure 1. Proton NMR Spectrum Expanded Region (4.8-0.2 ppm) for a Vehicle Control Sample that Contains Media on which Cells Were Grown. This region contains the peaks of importance for the chemometrics analysis described.




Figure 2. PC2 versus PC1 Scores Plot for PCA Involving 76 Spectra.
Label Key: RC = Ricin Cells samples, RO = Ricin Blanks, BK = Media
Blanks, VC = Vehicle Controls. Color Key: Batch 1 = blue,
Batch 2 = red, Batch 3 = green, Batch 4 = magenta.




Figure 3. PC3 versus PC4 Scores Plot for PCA Involving 76 Spectra. Key: RC = Ricin Cells samples, RO = Ricin Blanks, BK = Media Blanks, VC = Vehicle Controls.




Figure 4. PC2 versus PC1 Scores Plot for PCA Involving Batch 1 Spectra. Key: RC = Ricin Cells samples, RO = Ricin Blanks, BK = Media Blanks, VC = Vehicle Controls.




Figure 5. PC2 versus PC1 Scores Plot for PCA Involving Batch 4 Spectra.
The numbers correspond to the ricin concentrations used for the
Ricin Cells samples, where #1 is the lowest concentration and VC is the
vehicle control sample. The ellipses are at the 0.95 confidence level.